Composition in %

Au	Pt	Pd	Rh	Ir	Ru	Ag	Cu	Ga	Fe	In	Sn	Zn
		38,8			0,1	55,0				1,1	2,0	3,0

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Directions for Use

0 Range of applications

This alloy, as a conventional casting and PressOver alloy, is suitable for the stated range of applications.

The alloy Simidur PressX can additionally be fired as an all-purpose alloy with low-fusing, high-expanding veneering ceramics.

1 Waxing-up

Minimum wall thickness 0.3 mm for single crowns and min. 0.4 mm for bridge abutment crowns. Cut back frameworks anatomically to allow for the planned veneer and ensure stable shaping for larger span widths; avoid sharp transitions. Create large-area surfaces for planned soldering.

2 Spruing

For single crowns and bridges up to max. 3 units, direct spruing (min. 3 mm ø) or pear-shaped sprues are recommended. Generally speaking, the cast bar (approx. 4 – 5 mm ø) is fitted with suitable guide channels (3 – 4 mm ø) and connecting channels (approx. 3 mm x approx. 3 mm ø) and, if appropriate, cooling ribs (1 mm ø).

3 Investing

Line casting flask with ceramic fleece. Mix the recommended investment material (see data chart (A)) under vacuum and invest (follow manufacturer's instructions carefully).

4 Preheating

With conventional graduated preheating methods, hold first preheating step at 280 °C for 30–60 min; further preheating according to the instructions of the investment material manu-facturer. Once the final temperature is reached (see data table (B)), the holding time depending on the size of the flask is: 20 min (1x), 40 min (3x), 50 min (6x), 60 min (9x).

With speed preheating methods it is essential to follow the in-structions of the investment material manufacturer precisely.

5 Melting/Casting

See data table (C, D, E) for casting crucible, casting temperature and post-melting time, depending on the method used. Using a separate crucible for each alloy is recommended. With regard to re-using sprues, it is advisable only to use thoroughly cleaned casting residues and add at least half new material.

6 Finishing/Cleaning

After slow cooling until warm to the touch, devest and, if appropriate, blast with glass beads (approx. 50 µm) or dip (Wilavest Sonic Liquid). Finish framework with car-bide cutters. Then blast surfaces with pure aluminum oxide (approx. 110 µm) under a pressure of about 2 bar. Brush off the framework thoroughly under running water or blast with super-heated steam and degrease with ethyl alcohol.

7 Oxidizing

Oxidize without a vacuum; set temperature and time according to data table (F). Then repeat blasting and cleaning of the parts being overpressed.

8 Firing on PressOver opaquer

Application of the PressOver opaquer in two firings is recom-mended. The PressOver opaquer must be applied according to the instructions in the relevant directions for use. It is important to ensure adequate support of the metal framework.

9 Overpressing

When overpressing the metal framework fired with PressOver opaquer, the rules for fabricating pressed ceramic frameworks must be strictly adhered to: waxing up, minimum wall thickness, spruing, pressing parameters.

10 Joining

For soldering before overpressing or of non-overpressed cast frameworks, use solder as stipulated in the data table (H); soldering gap 0.05 to 0.2 mm. Allow soldered object to cool slowly after soldering.

For laser welding before overpressing, use welding wire as stipulated in the data table (J). Caution: pay attention to weld seam geometry.

11 Polishing of non-overpressed frameworks

Remove residual oxides and any flux residues by blast-polishing. Polish to a semi-gloss shine with rubber polishers (Polirol U), polish to a high-glaze finish with suitable polishing paste, then clean using ultrasonic cleaning equipment or by careful steam blasting.

Caution: Metal vapours and metal dust are harmful if inhaled; so always use ventilation and/or suitable protec-tive masks!

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